



**HAMDAM TEB
NOVIN LASER**

MEDICAL ENGINEERING



WINOLLA

Non-Contact Laser



- Four wavelengths (1064-940-808-755nm)
- Suitable for all skin types
- No skin contact
- Painless experience
- Detachable personal distance gauge
- Total power 3200W
- Optical power 1600W
- Large spot size (5.2×1cm)
- Internal chiller compressor + Dynamic air cooling
- Continuous operation and long lifespan
- 15.5 inch touch screen display

WINOLLA NON-CONTACT HAIR REMOVAL LASER INTRODUCTION



WEIGHTLESS COMFORT!

As a result of novel design of the handpiece holder, you won't experience any weight or fatigue during operation.



NON CONTACT LASER

Due to detachable personal distance gauge and elimination of gel the laser operation is carried out with complete hygiene, without transmission of any diseases. Additionally, with the use of dynamic air cooling, a pain-free and burn-free laser experience is ensured.

Winolla hair removal laser due to its four-wavelengths technology delivers the best results in the shortest possible time. Not only does the NIR wavelength help brighten the skin, but it also attacks blood vessels that nourish the hair and the capillaries.



Four-Wavelengths Energy Source



Detachable Personal Distance Gauge



Powerful Dynamic Air Cooling -20°C



Internal Chiller Compressor





OPTI MAX PRO

Non-Contact Laser



- Three wavelengths (1064-808-755nm)
- Suitable for all skin types
- No skin contact
- Painless experience
- Detachable personal distance gauge
- Total power 2800W
- Optical power 1200W
- Large spot size (4.8×1 cm)
- Internal chiller compressor + Dynamic air cooling
- Continuous operation and long lifespan
- 15.5 inch touch screen display

OPTIMAX PRO NON-CONTACT HAIR REMOVAL LASER INTRODUCTION

WEIGHTLESS COMFORT!

As a result of novel design of the handpiece holder, you won't experience any weight or fatigue during operation.



NON CONTACT LASER

Due to detachable personal distance gauge and elimination of gel the laser operation is carried out with complete hygiene, without transmission of any diseases. Additionally, with the use of dynamic air cooling, a pain-free and burn-free laser experience is ensured.

The Optimax Pro hair removal laser due to its three-wavelengths technology delivers the best results in the shortest possible time. In contrast to general Optimax, this product has a larger spot size and upgraded energy source, resulting in shorter operation time.



Combined Multi-Length Energy Source



Detachable Personal Distance Gauge



Powerful Dynamic Air Cooling -20c



Internal Chiller Compressor





OPTIMAX

Non-Contact Laser

2800W
1000W



- Three wavelengths (1064-808-755nm)
- Suitable for all skin types
- No skin contact
- Painless experience
- Detachable personal distance gauge
- Total power 2800W
- Optical power 1000W
- Large spot size (4×1cm)
- Internal chiller compressor + Dynamic air cooling
- Continuous operation and long lifespan
- 15.5 inch touch screen display



OPTIMAX NON-CONTACT HAIR REMOVAL LASER INTRODUCTION

Weightless Comfort!

As a result of novel design of the handpiece holder, you won't experience any weight or fatigue during operation.



NON CONTACT LASER

Due to detachable personal distance gauge and elimination of gel the laser operation is carried out with complete hygiene, without transmission of any diseases. Additionally, with the use of dynamic air cooling, a pain-free and burn-free laser experience is ensured.

Optimax hair removal laser due to its three-wavelengths technology delivers the best results in the shortest possible time.



Combined Multi-Length Energy Source



Detachable Personal Distance Gauge



Powerful Dynamic Air Cooling -20c



Internal Chiller Compressor





PARADISE **PRO**

Diode Laser Machine

1600



- Four wavelengths (1064-940-808-755nm)
- Suitable for all skin types
- Painless experience
- Optical power 1600W
- Equipped with sapphire crystal
- Large spot size (3.4×1.5 cm)
- Internal supercooling compressor
- Continuous operation and long lifespan
- 15.5 inch touch screen display



PARADISE PRO
HAIR REMOVAL LASER INTRODUCTION

WAVELENGTHS

IV

808 - 755
1064 - 940
(nm)



I

808 - 755
1064 - 940
(nm)



V

808 - 755
1064 - 940
(nm)



II

808 - 755
1064 - 940
(nm)



VI

808 - 755
1064 - 940
(nm)



III

808 - 755
1064 - 940
(nm)




4
WAVELENGTHS
TECHNOLOGY

Four-wavelengths technology (**755 Alexandrite** / **808 Diode** / **940 NIR** / **1064 Nd-YAG**) simultaneously radiates onto the skin and it delivers remarkable effectiveness on all skin and hair types in the fewest sessions.

The presence of the NIR wavelength helps brighten the skin, also it attacks blood vessels that nourish the hair and the capillaries.



Four-Wavelengths
1600W Source



Sapphire Crystal



Dual TEC Cooling



Internal Cooling
Compressor





PARADISE PLUS

Hair Removal Laser



- Three wavelengths (1064-808-755nm)
- Suitable for all skin types
- Painless experience
- Optical power 1200W
- Equipped with sapphire crystal
- Large spot size (3×1.5 cm)
- Internal supercooling compressor
- Continuous operation and long lifespan
- 15.5 inch touch screen display



PARADISE PLUS

HAIR REMOVAL LASER INTRODUCTION

VI

1064 - 808 - 755 (nm)



V

1064 - 808 - 755 (nm)



IV

1064 - 808 - 755 (nm)



III

1064 - 808 - 755 (nm)



II

1064 - 808 - 755 (nm)



I

1064 - 808 - 755 (nm)




DIODE TECHNOLOGY

Diode technology is one of the most popular and fastest method for hair removal. This Technology simultaneously radiates three wavelengths (Nd-YAG 1064 / Diode 808 / Alexandrite 755) onto the skin and attacks the melanin; additionally due to deep penetration, high precision and internal supercooling compressor, pain-free laser experience in shortest possible time delivers best results. This product also has a better energy source and larger spot size compared to Paradise Normal. These features will result in shorter operation time.



Three-Wavelengths
1200W Source



Sapphire Crystal



Dual TEC Cooling



Internal Cooling
Compressor





PARADISE **NORMAL**

Hair Removal Laser



- Three wavelengths (1064-808-755nm)
- Suitable for all skin types
- Painless experience
- Optical power 1000W
- Equipped with sapphire crystal
- Large spot size (2.5×1.5 cm)
- Internal supercooling compressor
- Continuous operation and long lifespan
- 15.5 inch touch screen display



PARADISE NORMAL HAIR REMOVAL LASER INTRODUCTION

VI

1064 - 808 - 755 (nm)



V

1064 - 808 - 755 (nm)



IV

1064 - 808 - 755 (nm)



III

1064 - 808 - 755 (nm)



II

1064 - 808 - 755 (nm)



I

1064 - 808 - 755 (nm)




DIODE TECHNOLOGY

Diode technology is one of the most popular and fastest method for hair removal. This Technology simultaneously radiates three wavelengths (**Nd-YAG 1064** / **Diode 808** / **Alexandrite 755**) onto the skin and attacks the melanin; additionally due to deep penetration, high precision and internal supercooling compressor, pain-free laser experience in shortest possible time delivers best results.



Three-Wavelengths
1000W Source



Sapphire Crystal



Dual TEC Cooling



Internal Cooling
Compressor






HIFU DEEP

High Intensity Focused Ultrasound



- Vaginal Handpiece
- Face / Body Handpiece
- 3 facial cartridges
- 2 body cartridges
- 2 vaginal cartridges
- Treats sagging face, wrinkles, fine lines
- Penetrates deep, middle and superficial layers
- Instant result
- Treats urinary incontinence
- Treats vaginal dryness
- No convalescence period
- 15.5 inch touch screen display

HIFU DEEP INTRODUCTION



HIFU

In this technology, focused ultrasonic waves target the dermis, damaged collagen fibers and new collagen fibers are produced through cellular stimulation, leading to the reduction of wrinkles and smoother skin.

SPEED OF ONE LINE 02 PER SECOND

3D Face & Body Cartridges

SPEED OF ONE LINE 01 PER SECOND

2D Vaginal Cartridges





TANINO EMS

Electromagnetic Muscle Stimulation



- 2 Flat handpieces
- 2 curved handpieces
- Simultaneous activity of all handpieces
- 30% Fat reduction
- 70% muscle growth
- Body contouring
- Chronic muscle and joint pain reduction
- Suitable for all ages
- No need for intense or prolonged exercise
- Blood circulation enhancement
- 30min of activity = 30000 sit-up
- 15.5 inch touch screen display



TANINO MUSCLE BUILDING MACHINE INTRODUCTION


Biceps sculpting and muscle growth
(Flat & Curved Handpiece)

Abdominal sculpting and muscle growth
(Flat Handpiece)

Flanks sculpting and muscle growth
(Curved Handpiece)

Lateral thigh sculpting and muscle growth
(Curved Handpiece)

Anterior thigh sculpting and muscle growth
(Curved Handpiece)




Triceps sculpting and muscle growth
(Flat Handpiece)

Waist sculpting and muscle growth
(Flat & Curved Handpiece)

Glutes sculpting and muscle growth
(Curved Handpiece)

Hamstrings sculpting and muscle growth
(Curved Handpiece)

Gastrocnemius sculpting and muscle growth
(Flat & Curved Handpiece)



EMS

In EMS, the focused electromagnetic waves instantly cause tension and extension of muscles, leading to muscle growth without intense exercise in the targeted area.

2 CURVED HANDPIECES

For curved areas

4 WARMING HANDPIECES

Blood circulation enhancement

2 FLAT HANDPIECES

For flat areas





ULTIMO PLUS

Slimming Machine




- 4 Cryo 360° body handpieces
- 1 Cryo 360° double chin handpiece
- 1 (40KHz) Cavitation handpiece
- 1 RF Body handpiece
- 1 RF Face handpiece
- 360° Freezing
- Simultaneous cryo activity
- First-session size reduction
- Internal supercooling compressor
- Temperature reduction down to -10°C
- 4.5 inch handpiece touch screen display
- 15.5 inch touch screen display



ULTIMO PLUS

SLIMMING MACHINE INTRODUCTION



CRYO
CAVITATION
RF

Cryolipolysis: Freezing of interstitial fluid causes fat cell death.

Cavitation: Ultrasonic waves break resistant fat.

RF: Radiofrequency waves eliminate fat located in the deeper layers of the body.



Internal Cooling
Compressor



RF Face



40KHz cavitation



RF Body





ULTIMO NORMAL

Slimming Machine




- 2 Cryo 360° body handpieces
- 1 Cryo 360° double chin handpiece
- 1 (40KHz) Cavitation handpiece
- 1 RF Body handpiece
- 1 RF Face handpiece
- 360° 360° Freezing
- 2 Simultaneous cryo activity
- First-session size reduction
- Internal supercooling compressor
- Temperature reduction down to -10°C
- 4.5 inch handpiece touch screen display
- 15.5 inch touch screen display



ULTIMO NORMAL

SLIMMING MACHINE INTRODUCTION



CRYO
CAVITATION
RF

Cryolipolysis: Freezing of interstitial fluid causes fat cell death.

Cavitation: Ultrasonic waves break resistant fat.

RF: Radiofrequency waves eliminate fat located in the deeper layers of the body.



Internal Cooling Compressor



RF Face



40KHz cavitation



RF Body





PANDORA **PLUS** LPG

Slimming Machine




- 1 LPG handpiece
- 1 (80KHz) cavitation handpiece
- 1 (40KHz) cavitation handpiece
- 1 (3MHz) ultrasonic handpiece
- 1 RF Body handpiece
- 1 RF Face handpiece
- First-session size reduction
- Face and Body lifting
- 10 inch touch screen display



PANDORA PLUS LPG

SLIMMING MACHINE INTRODUCTION



LPG

This feature is designed to combine both massage and suction, gently massaging the skin while helping to break down the accumulated fat. The roller's movement increases blood circulation in the area, which helps to break down fat cells and facilitate their elimination from the body.



40 / 80KHz
Cavitation




3 MHz
Ultrasonic



RF Body



RF Face





PANDORA **PLUS3C**

Slimming Machine




- 1 RF-Vacuum handpiece
- 1 (80KHz) cavitation handpiece
- 1 (40KHz) cavitation handpiece
- 1 (3MHz) ultrasonic handpiece
- 1 RF Body handpiece
- 1 RF Face handpiece
- First-session size reduction
- Face and Body lifting
- 10 inch touch screen display



PANDORA PLUS 3C

SLIMMING MACHINE INTRODUCTION



Pandora Plus 3C utilizes not only body and face RF but also incorporates two cavitations and one ultrasonic with various frequencies.

40KHz cavitation : suitable for resistant fat in deep layers

80KHz cavitation : suitable for resistant fat in superficial layers

3 MHz ultrasonic : suitable for lifting sensitive areas like face and neck



40 / 80KHz
Cavitation



3 MHz
Ultrasonic



RF Body



RF Face



BEAUTY DEVICES



LASER, SLIMMING AND REJUVENATION

EMAIL@HAMDAMTEB.COM
0912 9731884
021 26756724
WWW.HAMDAMTEB.COM